




COEZYME™

DUAL ENZYMATIC DETERGENT CONCENTRATE

PRODUCT DESCRIPTION

COEZYME is a highly concentrated, dual enzymatic detergent. Available in self-metering quart and half gallon pump bottles, COEZYME is easily mixed with tap water. Used as directed, COEZYME prevents solidification of blood and bioburden, and disintegrates these protein substances (i.e. blood, tissue and mucous) if they have dried.



BENEFIT STATEMENT

Immersion in COEZYME minimizes manual handling/scrubbing of contaminated items, reducing the risk of percutaneous injury. IONIC SURFACTANTS help the powerful solution to access and clean hard-to-reach areas. The synergistic enzymes dissolve and lift proteins, and the low-sudsing, neutral pH detergent component removes the dissolved particulates.

ONE OUNCE MAKES ONE GALLON!

PRODUCT CODES

FEATURES

- Detergent has neutral pH
- Low sudsing
- Powerful enzymatic action that dissolves and lifts protein
- Dual purpose
- Pre-calibrated pump
- Highly concentrated

BENEFITS

- Guards against instrument deterioration
- Compatible for use in ultrasonic cleaner
- If using as a presoak, manual cleaning is minimized and made easier
- If used in an ultrasonic cleaner, manual cleaning is virtually eliminated
- Staff exposure to potentially infectious contaminants via manual cleaning is kept to a minimum
- Multiplicity/simplicity of use
- Easy and accurate mixing
- Saves storage space and shipping costs



GC AMERICA INC.
3737 WEST 127TH STREET
CHICAGO, ILLINOIS 60658
TEL (708) 597-0900
FAX (708) 371-5103