



GC AMERICA INC.
3737 WEST 127TH STREET
ALSIP, ILLINOIS 60803
TEL (708) 597-0900
FAX (708) 371-5103

Contact:
Marketing Communications Specialist
GC America Inc.
3737 W. 127th Street
Alsip, IL 60803
708-897-4051

FOR IMMEDIATE RELEASE

GC America Introduces: GC Fuji TEMP LT™

Alsip, IL – February 2011 - GC America is very proud to introduce GC Fuji TEMP LT, the first of its kind conventional Glass Ionomer Temporary Cement that will be available in March.

GC Fuji TEMP LT is indicated for:

- Temporary (provisional) cementation of all types of all-ceramic, resin and metal-based crowns and bridges, as well as implant abutments.

GC Fuji TEMP LT provides the following features and benefits:

- Excellent handling with no-runny consistency
- Easy to remove excess cement which helps prevent remnants which can cause peri-implantitis
- Low film thickness (6 microns) (CAD CAM made prosthetics are very precise fitting)
- Radiopacity to make overhangs visible (290%)
- High bond strength and mechanical properties for cases where high retention is needed
- Excellent thermal insulation (prevents sensitivity)
- Excellent biocompatibility
- High Fluoride Release and recharge
- No adverse effect on final cementation which is common with eugenol cements
- Easy clean-up once you remove the provisional to lute the final appliance

For more information on GC America's products, please call Customer Service at 800-323-7063 or visit www.gcamerica.com

###

