



G-CEM LinkAce™

SELF-ADHESIVE RESIN CEMENT




ONE STEP FOR A HIGH BOND

TO PRETREATED

- PFM
- ZIRCONIA
- LITHIUM DISILICATE



Tensile Bond Strength to Dentin¹



TOUCH of GENIUS™

- **Optimal Self-cure Mode:** Innovative initiator system (chemical)
- **High Bond Durability to Pretreated Zirconia in ONE Step:** Proprietary phosphate monomers
- **Exceptional Color Stability:** Very low water sorption and HEMA-free
- **Unsurpassed Wear Resistance:** Small homogenously distributed particles
- **Easy Excess Cement Removal:** Only 1-2 second tack cure

A Winning Hand for all Indirect Restorations







1 HIGH BOND IN ONE STEP





Tensile Bond Strength to Pretreated Zirconia[†]



Tensile Bond Strength to Pretreated e.max[†] (Lithium disilicate)[†]



Tensile Bond Strength to Dentin / Self-cure Mode Only[†]



†Not a registered trademark of GC America Inc. 1GCC Internal data on file.



G-CEM LinkAce™ Packages
 Contains: 2 syringes (4.6 g - 2.7 mL each), 15 Automix Tips Regular, 5 Automix Tips for Endo.
 004857 AZ 004858 Translucent 004859 AO3
 004860 BO1

G-CEM LinkAce™ Accessories
 002902 20 GC Automix Tips (Regular) 002903 10 GC Automix Tips (Endo)