



THE IDEAL CONVERGENCE OF ADAPTATION, STRENGTH & ESTHETICS

Gænial™ BULK Injectable

Injectable High-Strength Ultra-Fine Particle Composite

Highly shapeable and easy to place, yet stays where you put it!

- Easy one-step filling – up to 4mm depth of cure with no capping or veneering required
- Full coverage silane coating of ultra-fine particles facilitates excellent adaptation to the cavity preparation
- Provides extraordinary wear-resistance and gloss retention
- Low volumetric shrinkage and low shrinkage stress helps ensure optimal marginal integrity







A restorative composite that is ideal for bulk filling to the occlusal surface of large cavities, by virtue of its injectable viscosity and excellent physical properties even superior to paste type composites*.

- Injectable composite with excellent viscosity and thixotropy for ideal bulk filling
- Can achieve highly durable restoration via high wear resistance and flexural strength (unlike competitive composite)
- Excellent polishability and gloss retention with self-polishing ability
- Guaranteed depth of cure of 4mm which allows for bulk filling in one layer

Class	Cavity Depth	Product Recommendation
I / II	Up to 2 mm	G-ænial™ Flo X or
		G-ænial™ Universal Injectable
I / II	Up to 4 mm	G-ænial™ BULK Injectable
I / II	4-6 mm	G-ænial™ BULK Injectable & G-ænial Sculpt®
I / II	6 mm +	G-ænial™ BULK Injectable & G-ænial Sculpt®

G-ænial™ BULK Injectable provides the right combination of strength and esthetics to provide a durable restoration that looks great at every recall!

- High flexural strength to protect chipping or fracture
- Outstanding wear resistance to provide ultimate durability, with excellent color and gloss retention
- Low volumetric shrinkage and low shrinkage stress to preserve marginal integrity and protect against recurrent caries



Courtesy of Jeff Lineberry, DDS, Mooresville, NC

G-ænial™ BULK Injectable Refills

Contains: One 3.4g (2.0mL) syringe, 20 Dispensing Tips III (plastid) and one light protective cap.
012394 A1 012395 A2

G-ænial Sculpt® Syringe Refills

Contains: One 2.0mL Syringe.
009152 A1 009153 A2 009154 A3 009155 A3.5 009156 A4
009168 AE 009161 AO1 009162 AO2 009163 AO3 009157 B1
009158 B2 009159 B3 009164 BW 009160 C3 009166 CV
009167 JE 009165 XBW

G-ænial Sculpt® 20-Unitip Refills

Contains: 20 0.16 mL unitips.
009169 A1 009170 A2 009171 A3 009172 A3.5 009173 A4
009174 B1 009175 B2 009176 B3 009177 C3

G-ænial Sculpt® 10-Unitip Refills

Contains: 10 0.16mL unitips.
009185 AE 009178 AO1 009179 AO2 009180 AO3 009180 BW
009183 CV 009184 JE 009182 XBW