



BRING THE X FACTOR TO POSTERIOR RESTORATIONS

G-aenial™ Flo X

Light-Cured Radiopaque Flowable Composite

- Homogenously dispersed ultra-fine barium-glass particles provide excellent handling
- Optimal viscosity for lining; flows easily and ensures excellent adaptation to cavity walls
- High radiopacity, easily detectable on radiograph even at low thickness







Featuring high radiopacity, easily detectable on radiograph even at low thickness, excellent handling properties, and optimal durability, G-ænial™ Flo X is the perfect partner to G-ænial™ Universal Injectable.

- Indicated for use in small restorations and as a cavity liner
- Highly radiopaque and highly flowable
- Provides easy placement and flow into preparations for excellent handling

Class	Cavity Depth	Product Recommendation
I / II	Up to 2 mm	G-ænial™ Flo X or
		G-ænial™ Universal Injectable
I / II	Up to 4 mm	G-ænial™ BULK Injectable
I / II	4-6 mm	G-ænial™ BULK Injectable & G-ænial Sculpt®
I / II	6 mm +	G-ænial™ BULK Injectable & G-ænial Sculpt®

G-ænial™ Flo X is a highly versatile flowable composite, with high flexural strength and excellent wear resistance.



G-ænial™ Flo X Refills
 Contains: One 3.6g (2.0mL) syringe, 20 Dispensing Tips III (plastic) and one light protective cap.
 008317 A1
 008318 A2
 008319 A3
 008320 A3.5
 008321 A4
 008322 AO2
 008323 AO3
 008324 CV

