



Contact:

Brendan Steidle
Product Marketing Manager
3737 W. 127th Street
Alsip, IL 60803
708-926-3007
brendan.steidle@gc.dental

GC America Introduces:

G-aenial™ Flo X

Light-Cured Radiopaque Flowable Composite

Alsip, IL – August 2018 – GC America Inc. is pleased to announce the **October 2018** launch of G-aenial™ Flo X composite.

G-aenial™ Flo X is a low-viscosity flowable composite which is perfect for all your posterior restorations! It has been designed to ensure excellent sealing and adaptation to cavity walls to effectively avoid porosities, preserve marginal integrity, reach even the most narrow and difficult to access cavities, and offer a minimally invasive treatment of caries to the patient. It offers a great balance between strength and viscosity in a flowable, with very high wear-resistance. G-aenial™ Flo X features a higher radiopacity, even at low thickness, to ensure perfect x-ray visibility and easy diagnostic of secondary caries for the dentist. G-aenial™ Flo X is mainly indicated as a liner or base, but may also be used for filling undercuts, tunnel restorations, fissure sealing, and Class V restorations.

G-aenial™ Flo X is available in a 1 x 3.6g (2.0 mL) syringe, in A1, A2, A3, A3.5, A4, AO2, AO3, and CV shades.

For more information on GC America and its complete product line, please visit our corporate website at www.gcamerica.com.

#####