

**GC QUALITY**  
**+ Strontium Glass + FLUORIDE**  
**+ Tack Cure + SYRINGE**  

---

**= Evolve**



# THE SOLUTION IS EVOLUTION

Introducing GC FujiCEM® Evolve in New Syringe Delivery System

**GC FujiCEM® Evolve**  
Resin Modified Glass Ionomer Cement

Building on the legacy GC FujiCEM® 2, new GC FujiCEM® Evolve takes the luting process to the next level!

- Superior bond strength to zirconia than a leading competitor<sup>1</sup>
- New tack-cure feature reduces cleanup time to seconds
- High radiopacity for easy visualization
- Ideal for cementation of zirconia, PFM and lithium disilicate restorations
- Rechargeable fluoride release and moisture tolerance

Try new GC FujiCEM® Evolve and join the evolution.



<sup>1</sup> Fujimoto et al, Bonding Durability of Novel Resin-Modified Glass-Ionomer Luting Cement to Zirconia, J Dent Res 98 (Spec Iss A): #1315, 2019.



## Tack cure (3 seconds) for easy excess removal



Courtesy of Dr. Javier Tapia Guacix (Spain)

## Glass ionomer advantages



Moisture tolerant – no need for rubber dam



Technique tolerant with no need for primer

Courtesy of Dr. Jeff Lineberry (USA)

## GC FujiCEM® Evolve is especially indicated for:



High strength ceramic (zirconia and lithium disilicate) crowns and bridges.



Metal and porcelain fused to metal (PFM) crowns and bridges.




Tack-cure ability with good working time (2'15) and fast setting (4'30)



Virtually no post-operative sensitivity


Courtesy of Dr. Jeff Lineberry (USA)

### Superior Bond Strength to Zirconia Than Many Competitors (MPa)




Source: GCC R&D

### Superior Radiopacity Than Many Competitors (AI%)




Source: GCC R&D

## How to lock?



## How to unlock?



### 012949 GC FujiCEM® Evolve Triple Pack Handmix

Contains: Three 9.2g / 5.0mL syringes and one mixing pad.

### 012950 GC FujiCEM® Evolve Triple Pack Automix

Contains: Three 9.2g / 5.0mL syringes and 45 GC Push and Click Tips (Regular).

### 012948 GC FujiCEM® Evolve Single Pack

Contains: One 9.2g / 5.0mL syringe and 15 GC Push and Click Tips (Regular).

### 012954 GC Push and Click Tips (Regular x 15)

### 012955 GC Push and Click Tips (Endo x 15)

