



# NOGENOL™

## NON-EUGENOL TEMPORARY CEMENT

### PRODUCT DESCRIPTION



**NOGENOL** is indicated for temporary cementation of plastic crowns or other provisional restorations. **NOGENOL** possesses adequate adhesive strength to resist dislocation and leakage, yet allows for easy removal of crowns with KY Pliers. Ideal for use with ALIKE™ or UNIFAST LC™ or UNIFAST LC™ Temporary C & B Acrylic.

### FEATURES AND BENEFITS

<b>Contains no eugenol</b>	• Does not burn soft tissue
<b>Low solubility</b>	• Provides excellent marginal seal
<b>Will not craze plastic / polycarbonate crowns</b>	• Extended longevity of temporary use
<b>Will NOT affect acid-etch procedures</b>	• Can be used for temporary inlays / onlays
<b>EIGHTEEN MONTH shelf life</b>	• Will not spoil in slow-average rate of use

### PRODUCT CODES

**136402 NOGENOL Temporary Cement, Standard Package**  
 Contains: 2 Tubes Base (30 g ea.), 1 Tube Catalyst (65 g)



**GC America Inc.**

*Advancing the Art and Science of Dentistry*

© 1998 GC America Inc., 3737 West 127th St., Alsip, Illinois 60803  
 Printed in U.S.A.