



# OMNIFLEX™

## POLYSULFIDE IMPRESSION MATERIAL

### PRODUCT DESCRIPTION

OMNIFLEX is a SINGLE MIX, all-purpose rubber base impression material. Recommended for all types of crowns, bridges, inlays and temporary denture relines. In more than 20 years of clinical experience, OMNIFLEX has consistently provided precise, detailed impressions. Thirty-six month shelf life.



### FEATURES AND BENEFITS

Unique consistency is suitable for syringe and tray - single mix

- Saves time, simplifies ordering.

Patented copper catalyst system

- "Snap" set in the mouth. Does not stain.

Excellent tear strength

- Resists tearing in deep undercuts.

High elasticity

- Easy to remove from stone cast.

Blue-green color, mint flavor

- Enhanced patient acceptance.

Available in 2 set times:

Regular - 3½ minute work time  
4 minute set time

Fast Set - 2½ minute work time  
3 minute set time

- Suitable for a wide range of impression techniques.

### PRODUCT CODES

	<b>Standard Package</b>			<b>OMNIFLEX—Catalyst Only (45 ml)</b>
	Contains: 1 Base (90 ml), 1 Catalyst (45 ml)	<b>137011</b>		<b>Regular Set</b>
	Adhesive (15 ml)	<b>137311</b>		<b>Fast Set</b>
<b>137001</b>	<b>Regular Set</b>			
<b>137301</b>	<b>Fast Set</b>			
	<b>Bulk Package</b>			<b>OMNIFLEX—Base Only (90 ml)</b>
	Contains: 36 Base (90 ml ea.),	<b>137021</b>		<b>Regular Set</b>
	36 Catalyst (45 ml ea.)	<b>137321</b>		<b>Fast Set</b>
<b>137036</b>	<b>Regular Set</b>	<b>133901</b>		<b>COE Adhesive Package for all polysulfides</b>
<b>137336</b>	<b>Fast Set</b>			(2 bottles COE Adhesive, 1 bottle thinner)
		<b>133911</b>		<b>COE Adhesive 15 ml (½ oz.)</b>
		<b>133912</b>		<b>COE Adhesive 15 ml (½ oz.) Package of 12</b>
		<b>133921</b>		<b>COE Adhesive Thinner 15 ml (½ oz.)</b>



**GC America Inc.**

*Advancing the Art and Science of Dentistry*

© 1997 GC America Inc., 3737 West 127th St., Alsip, Illinois 60803  
Printed in U.S.A.