


EQUIA Forte® HT

A glass hybrid restorative system for long-term restorations.

MI Paste® ONE

Anti-cavity toothpaste with bio-available calcium, phosphate and fluoride.

RESTORE.
RECHARGE.
REMINERALIZE.



TWO FOR WON!

- High fluoride release at tooth-restorative interface by EQUIA Forte® HT
- Optimal marginal seal for long-term resistance to microleakage by EQUIA Forte® HT
- Calcium and fluoride rechargeable EQUIA Forte® HT restorations
- Remineralization and sensitivity reduction by MI Paste® ONE
- Exceptional treatment for pediatric, elderly and high-carries risk patients








Don't Just Restore - Recharge!

You may know that glass ionomers and glass hybrids offer a number of advantages. They work well in a moist environment, requiring no isolation. They form a chemical bond with dentin that provides excellent marginal seal. Plus, they offer a sustained release of fluoride over a long period, helping to heal and prevent caries lesions. But what you might not know is that when they are paired with the new MI Paste® ONE, you can extend the therapeutic remineralizing benefits to the whole mouth. Transform your restorative procedures into a remineralizing tour de force with EQUIA Forte® HT and new MI Paste® ONE.

In addition to fluoride, GC glass ionomer and hybrid products release strontium and calcium ions. Strontium ions aid in remineralization and help increase hardness of enamel and its acid resistance. A recent research study showed that GC Fuji TRIAGE® and EQUIA Forte® absorb and re-release calcium and fluoride from MI Paste® ONE.*

Therefore, these restorations can be recharged with calcium and fluoride when used in conjunction with MI Paste® ONE.



EQUIA Forte® HT Intro Kit

Contains: 48 EQUIA Forte® HT Capsules (Powder 0.4g, Liquid: 0.13g (0.10mL) per capsule), one EQUIA Forte® Coat (4mL) bottle, 1 CAVITY CONDITIONER (GC) (5.7mL).

453001 A1 453002 A2 453003 A3 453004 A3.5 453006 B1
453007 B2 453008 B3 453009 C4 453005 Assorted

EQUIA Forte® HT Fil Refill

Contains: 48 EQUIA Forte® HT Capsules (Powder 0.4g, Liquid: 0.13g (0.10mL) per capsule).

453010 A1 453011 A2 453012 A3 453013 A3.5 453014 B1
453015 B2 453016 B3 453017 C4 453000 Assorted

452500 EQUIA Forte® Coat Bottle Refill

Contains: One EQUIA Forte® Coat (4mL) bottle.

452505 EQUIA Forte® Coat Unit Dose Refill

Contains: 50 EQUIA Forte® Coat (0.1mL each) unit doses, 50 Micro-Tip Applicators and one Applicator Holder.

437000 MI Paste® ONE - Fresh Mint

Contains: 10 tubes (46g each).

000110 CAVITY CONDITIONER (GC)

Contains: One (5.7mL) bottle.

437555 GC CAPSULE APPLIER III

Contains: One capsule applicator.