



# Transform the way you work

Start injecting with our strongest  
direct restorative ever\*

**G-ænial™ Universal  
Injectable**



**GC**

An injectable high-strength ultra-fine particle composite with ideal viscosity, handling and adaptation characteristics that may be used for long-lasting esthetic restorations.

The homogeneously dispersed ultra-fine barium-glass particles formula (150nm particle size, 69% wt filler loading) provides high flexural strength and wear resistance for a durable and long lasting restoration.



Metal tip may be bent with fingers if necessary to better access posterior distal region.

Cavity Depth	Product Recommendation
Up to 2 mm	<b>G-ænial™ Flo X</b> or
	<b>G-ænial™ Universal Injectable</b>
Up to 4-6+ mm	<b>G-ænial™ BULK Injectable with G-ænial Sculpt®</b> (veneer layer optional)

G-ænial™ Universal Injectable features high flexural strength, for protection against chipping and fracture, and outstanding wear resistance for excellent durability and gloss retention.



#### G-ænial™ Universal Injectable Syringe Refills

Contains: Two 1.7g (1.0mL) syringes, 10 long dispensing tips and 10 short dispensing tips.

012364 A1	012365 A2	012366 A3	012367 A3.5
012368 A4	012369 B1	012370 B2	012371 CV
012372 CVD	012373 AO1	012374 AO2	012375 AO3
012376 JE	012377 AE	012378 XBW	012379 BW



**Now Available  
in Unitips!**

#### G-ænial™ Universal Injectable Unitip Refills

Contains: 15 0.16mL unitips.

012348 A1	012349 A2	012350 A3	012351 A3.5
012352 A4	012353 B1	012354 B2	012355 CV
012356 CVD	012357 AO1	012358 AO2	012359 AO3
012360 JE	012361 AE	012362 XBW	012363 BW

425000 Unitip **APPLIER II**