

A MODULAR SYSTEM

GC Initial Bleach Shades are part of the GC Initial modular ceramic system and are available as separate sets. Each GC Initial Bleach Set contains different components and by adding them together you can create increasingly aesthetic restorations.

For more detailed information about the different set combinations please request our brochure.



Contents of GC Initial Bleach Sets

	GC Initial MC GC Initial LF GC Initial Ti	GC Initial AL	GC Initial Zr
Opaque/Liner/Frame Modifier	Opaque Modifier OM-1	BL-Lin	BL-FM
Dentin 20g	BLD-1	BLD-1	BLD-1
	BLD-2	BLD-2	BLD-2
	BLD-3	BLD-3	BLD-3
Bleach Enamel 20g	BL-E	BL-E	BL-E
Enamel Opalescence 20g	EOP-1	EOP-1	EOP-1
Fluo Dentine 20g	FD-91	FD-91	FD-91
Modelling Liquid 2 x 25ml	✓	✓	✓
Opaque Liquid 1 x 25ml	✓		
Colour shade tabs of all shades	✓	✓	✓

GC EUROPE N.V.
Head Office
Interleuvenlaan 13
B - 3001 Leuven
Tel: +32.16.39.80.50
Fax: +32.16.40.02.14
E-mail: info@gceurope.com
www.gceurope.com

AVAILABLE SEPARATELY



All bleach colours are available in 20g refill bottles.

Recommended complementary powders from the GC Initial range

For coloured individualisation, Bleach Dentin, BL-E, EOP-1, FD-91 can be mixed with all powders from the respective GC Initial assortment, below are some recommendations. They are available in all five GC Initial ceramic systems as refills or as part of modular sets.

ODM1	white Opaque Dentin
EOP-2	milky white Opal Enamel
EOP-3	light bluish Opal Enamel
EOP-4	light grey Opal Enamel
TO	opalescent Translucent Ceramic
CLF	fluorescent Translucent Ceramic
EO-15/16	yellowish white Opaque Enamel
ST-30	white, Translucent Shoulder Material
SO-37	white, Opaque Shoulder Material

For more detailed information, please read the corresponding GC Initial technical manual

GC UNITED KINGDOM Ltd.
22-23, Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel: +44.1908.218.999
Fax: +44.1908.218.900
E-mail: info@uk.gceurope.com
www.uk.gceurope.com

GC
FIRST IS QUALITY

Z L L F GB 7 40 04/05

Create your world with *initial*

The all-round ceramic system
for every indication

new

Bleach Shades

GC Initial Bleach Set MC • GC Initial Bleach Set LF • GC Initial Bleach Set Ti
GC Initial Bleach Set AL • GC Initial Bleach Set Zr



The GC Initial System

The GC Initial System is a line of six specialised ceramics that can be used to create metal-ceramic and full-ceramic restorations for every indication using any fabrication process and framework (for more information on the complete GC Initial system please request our brochure). To extend the unique harmonised shading system of GC Initial and follow the trend for very light tooth colours, GC has introduced GC Initial Bleach Shades, a unique colour set that allows you to match even the lightest tooth shades whatever ceramic restoration system you use.

How to use GC Initial Bleach Shades

Natural enamel reacts individually to bleaching therapies so a translucency assessment should be made six weeks after the final bleaching to ensure colour stability. GC Initial Bleach Shades have been specifically developed so they can be used with the five GC Initial ceramic systems (see table on the right). It should be noted that when it comes to colour-matched restorations, white opaque/opaque dentin should only be applied when minimal layering thickness is available to respect the original tooth chroma.



www.gcinitial.gceurope.com

GC
FIRST IS QUALITY

GC Initial MC/LF/Ti	GC Initial AL	GC Initial Zr
TO MASK COPINGS	AS A BASE FOR COPINGS	
Powder Opaque OM-1 Masks metal copings and can be used pure or in combination with standard opaques	AL Liner BL-Lin A base for aluminium oxide frame materials. It has a high level of fluorescence and a high wettability for easy application	Zr Frame Modifier BL-FM A base for zirconium oxide frame materials. It has a high level of fluorescence and a high wettability for easy application
TO MATCH DENTIN		
Bleach Dentin Three shades that vary in their value/chroma level. Use for layering instead of standard dentin. GC Initial BLD - 1: bleach dentin light (more white than the Vita® Classical A1 dentin shade) GC Initial BLD - 2: bleach dentin white (more white than the Vita Classical B1 dentin shade) GC Initial BLD - 3: bleach dentin Xwhite, extreme white Fluorescent Opaqus Dentin GC Initial Fluo dentin-91 (FD-91): white, high fluorescent opaques dentin. To cover the refraction of light from substructures/layering build-up. Is appropriate for bright white Mamelon structures.		
TO MATCH ENAMEL		
Bleach Enamel GC Initial BL-E: standard enamel shade. Light white base with an average translucency level. GC Initial Enamel Opal-1 (EOP-1): a very bright, strongly opalescent enamel shade with a slight cold-white effect. Frequently occurring colour effect with repeated and/or very strong bleaching therapies.		

*Vita is a registered trademark of Vita Zahnfabrik, Bad Säckingen.