


ALIKE™

PRODUCT DESCRIPTION



ALIKE is a self curing acrylic resin for constructing temporary crowns, bridges and splints. **ALIKE** is available in 6 shades cross referenced to BIOFORM® and VITA® ranges. Stain resistant **ALIKE** trims easily and will not clog burs.



FEATURES AND BENEFITS

Available in 6 shades	• Matches a wide variety of patients, ages
Minimal shrinkage	• Assures accurate fit, snug margins
High strength (surface hardness: 24,750 psi)	• Resists fracture, shear stresses
5-6 minute snap set	• Easy and convenient to use
Low water sorption (0.45 mg / cm ²)	• Reduced discoloration

PRODUCT CODES

340501 ALIKE—Professional Package
Contains: 6 X 1 1/2 oz. (45 g) Bottles Powder
(Shades B-59, 62, 65, 67, 77, 81)
4 oz. (118 ml) Liquid, Liquid Dropper

ALIKE Powder Refills: 1 1/2 oz. (45 g)

- 340511** Shade #59 (VITA® B-1)
- 340512** Shade #62 (VITA® A-1)
- 340513** Shade #65 (VITA® A-2)
- 340514** Shade #67 (VITA® A-3)
- 340515** Shade #77 (VITA® B-3)
- 340516** Shade #81 (VITA® A-3.5)

ALIKE Powder Refills: 6 oz. (170 g)

- 340521** Shade #59
- 340522** Shade #62
- 340523** Shade #65
- 340524** Shade #67
- 340525** Shade #77
- 340526** Shade #81

340591 ALIKE Liquid Refill 4 oz. (118 ml)