

# EXA MATERIALS GUIDE



## For EXAMIX, EXAFLEX and EXAFAST Products Vinyl Polysiloxane Material Usage Guide

Technique	Tray Material	Wash or Syringe Material
<b>Single Mix</b> <i>also Monophase, Medium Viscosity, Single Phase or Single Impression</i>	<b>Monophase</b> +	<b>Monophase</b>
<b>Double Mix</b>	<b>Heavy Body</b> OR <b>Monophase</b> +	<b>Regular</b> OR <b>Injection</b>
<b>Double Arch</b> <i>-Three Function</i>	<b>Heavy Body</b> +	<b>Regular</b> OR <b>Injection</b>
	OR <b>Monophase</b> +	<b>Monophase</b>
<b>Putty-Wash</b> <i>2-Step</i>	<b>Putty</b> +	<b>Regular</b> OR <b>Injection</b>
<b>Putty-Wash</b> <i>1-Step also Dimple Technique</i>	<b>Putty</b> +	<b>Regular</b> OR <b>Injection</b>