

## G-CEM LinkForce®

(GC America)

The **G-CEM LinkForce®** system features strong adhesion to all dental substrates in all indirect restorations, including crowns and veneers.

The system is simple to use and contains all components needed for successful placements. In addition to superior adhesion, **G-CEM LinkForce** has excellent color stability and tooth-like fluorescence. It may be used in both light-cure and self-cure modes. **G-CEM LinkForce** may be used in conjunction with many restorative materials, including stacked feldspathic ceramic, pressed leucite ceramic, and lithium disilicate. The **G-CEM LinkForce** cement comes in four shades and is available with corresponding try-in pastes in each shade.

