

G-Premio BOND™

(GC America)

G-Premio BOND™

is a universal, 8th generation bonding agent that is compatible with total-etch, self-etch and selective-etch techniques, providing excellent versatility. It is perfectly adapted to all direct restorations and can also be used to repair indirect restorations without the use of primer. **G-Premio BOND** can also be used in combination with a silane when repairing glass or hybrid ceramic, and is also ideal for hypersensitivity. It is available in a bottle (5 mL) with an easy-to-use flip-top cap and a silicone cover to enable more precise drop dispensing, and a 50 unit dose pack (0.16 mL each).

