

GC FujiCEM 2 (GC America)



GC FujiCEM 2 is a second-generation resin-modified glass ionomer luting cement in paste delivery. It is indicated for a myriad of indirect restorations, including inlays, onlays, crowns and bridges, as well as materials made of metal, resin, lithium disilicate and zirconia. Powered by *F2 Flex Fuse Technology*, which incorporates high-elastic crosslinking monomers with a modified filler-surface treatment, it improves indirect restorations by offering increased strength, high fluoride release, low-film thickness and excellent marginal integrity.