

GC PLIERS


GC America
www.gcamerica.com



Clinician: Dr. Sabiha S. Bunek



Emery powder



Problem:

Many dentists and assistants use a pair of hemostats for removing provisional restorations. They don't grip crowns well and may easily chip or fracture restorations.

Solution:

GC PLIERS are well-constructed stainless steel forceps designed to easily remove and handle indirect restorations. The working ends are rasped, so they grab into restorations well. They also come with rubber pads that can be snapped on the working ends to provide cushion and grip when needed. The pliers effectiveness is enhanced when paired with an abrasive powder.