

Item ID	Item Number	Tray Name	WIDTH						LENGTH						HEIGHT					
			Target		Tolerance		Range		Target		Tolerance		Range		Target		Tolerance		Range	
			Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric	Inches	Metric
COE 1	260011	Perf. Regular Tray, #1	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 68 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17 - 27 mm
COE 4	260041	Perf. Regular Tray #4	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68 - 78 mm	2 6/16	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 7	260071	Perf. Regular Tray #7	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68 - 78 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	11/16	17 mm	±3/16	±5 mm	8/16 to 14/16	12-22 mm
COE 9	260091	Solid Regular Tray, #109	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57 - 67 mm	1 7/8	48 mm	±3/16	±5 mm	1 11/16 to 2 1/16	43-53 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 20	260201	Perf. Regular Tray #20	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 21	260211	Perf. Regular Tray #21	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 1/4	57 mm	±3/16	±5 mm	2 1/16 to 2 7/16	52-62 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 22	260221	Perf. Regular Tray #22	2 15/16	75 mm	±3/16	±5 mm	2 12/16 to 3 2/16	70 - 80 mm	2 1/8	54 mm	±3/16	±5 mm	1 15/16 to 2 5/16	49-59 mm	13/16	21 mm	±3/16	±5 mm	10/16 to 1	16-26 mm
COE 24	260241	Solid Regular Tray, #124	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	1 15/16	49 mm	±3/16	±5 mm	1 12/16 to 2 2/16	44-54 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 30	260301	Perf. Partial Tray #30	1 3/4	44 mm	±3/16	±5 mm	1 9/16 to 1 15/16	39 - 49 mm	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57-67 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 31	260311	Perf. Partial Tray #31	1 3/4	44 mm	±3/16	±5 mm	1 9/16 to 1 15/16	39 - 49 mm	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57-67 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 32	260321	Perf. Partial Tray #32	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	1 1/2	38 mm	±3/16	±5 mm	1 5/16 to 1 11/16	33-43 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 33	260331	Perf. Partial Tray #33	2 5/8	67 mm	±3/16	±5 mm	2 7/16 to 2 13/16	62 - 72 mm	1 7/16	37 mm	±3/16	±5 mm	1 4/16 to 1 10/16	32-42 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE 41	260411	Mcgowan Immediate Denture Tray #41	3 1/8	79 mm	±3/16	±5 mm	2 15/16 to 3 5/16	74 - 84 mm	2 5/8	67 mm	±3/16	±5 mm	2 7/16 to 2 13/16	62-72 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 42	260421	Mcgowan Immediate Denture Tray #42	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68 - 78 mm	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57-67 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 43	260431	Mcgowan Immediate Denture Tray #43	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68 - 78 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 44	260441	Mcgowan Immediate Denture Tray #44	2 11/16	68 mm	±3/16	±5 mm	2 8/16 to 2 14/16	63 - 73 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	13/16	21 mm	±3/16	±5 mm	10/16 to 1	16-26 mm
COE 45	260451	Mcgowan Immediate Denture Tray #45	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68 - 78 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 46	260461	Mcgowan Immediate Denture Tray #46	2 5/8	67 mm	±3/16	±5 mm	2 7/16 to 2 13/16	62 - 72 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE 47	260471	Mcgowan Immediate Denture Tray #47	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65 - 75 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	13/16	21 mm	±3/16	±5 mm	10/16 to 1	16-26 mm
COE 48	260481	Mcgowan Immediate Denture Tray #48	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	56 - 66 mm	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57-67 mm	13/16	21 mm	±3/16	±5 mm	10/16 to 1	16-26 mm
COE 80	260801	McGowan-Winkler Full Denture Tray #80	3 1/8	79 mm	±3/16	±5 mm	2 15/16 to 3 5/16	74 - 84 mm	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65-75 mm	9/16	14 mm	±3/16	±5 mm	6/16 to 12/16	9-19 mm
COE 81	260811	McGowan-Winkler Full Denture Tray #81	2 13/16	71 mm	±3/16	±5 mm	2 10/16 to 3	66 - 76 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	9/16	14 mm	±3/16	±5 mm	6/16 to 12/16	9-19 mm
COE 82	260821	McGowan-Winkler Full Denture Tray #82	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	9/16	14 mm	±3/16	±5 mm	6/16 to 12/16	9-19 mm
COE 84	260841	McGowan-Winkler Full Denture Tray #84	2 15/16	75 mm	±3/16	±5 mm	2 12/16 to 3 2/16	70 - 80 mm	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65-75 mm	1/4	6 mm	±3/16	±5 mm	1/16 to 7/16	1-11 mm
COE 85	260851	McGowan-Winkler Full Denture Tray #85	3 1/4	83 mm	±3/16	±5 mm	3 1/16 to 3 7/16	78 - 88 mm	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57-67 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE 86	260861	McGowan-Winkler Full Denture Tray #86	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	1/4	6 mm	±3/16	±5 mm	1/16 to 7/16	1-11 mm
COE 87	260871	Mcgowan-winkler Full Denture Tray #87	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 1/4	57 mm	±3/16	±5 mm	2 1/16 to 2 7/16	52-62 mm	1/4	6 mm	±3/16	±5 mm	1/16 to 7/16	1-11 mm
COE 88	260881	Mcgowan-winkler Full Denture Tray #88	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65 - 75 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE 89	260891	Mcgowan-winkler Full Denture Tray #89	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2 5/16	59 mm	±3/16	±5 mm	2 2/16 to 2 8/16	54-64 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE 99	260991	Perf. Partial Tray #99	1 1/4	32 mm	±3/16	±5 mm	1 1/16 to 1 7/16	27 - 37 mm	2 3/16	56 mm	±3/16	±5 mm	2 to 2 6/16	51-61 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE U1S	263471	Perf. STO-K Tray, U-1-S	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	56 - 66 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	7/16	11 mm	±3/16	±5 mm	4/16 to 10/16	6-16 mm
COE U2S	263472	Perf. STO-K Tray, U-2-S	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	3/4	19 mm	±3/16	±5 mm	9/16 to 15/16	14-24 mm
COE U3S	263473	Perf. STO-K Tray, U-3-S	3 1/4	83 mm	±3/16	±5 mm	3 1/16 to 3 7/16	78 - 88 mm	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65-75 mm	7/8	22 mm	±3/16	±5 mm	11/16 to 1 1/16	17-27 mm
COE L1S	263481	Perf. STO-K Tray, L-1-S	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57 - 67 mm	1 3/4	44 mm	±3/16	±5 mm	1 9/16 to 1 15/16	39-49 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE L2S	263482	Perf. STO-K Tray, L-2-S	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE L3S	263483	Perf. STO-K Tray, L-3-S	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57 - 67 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE L4S	263484	Perf. STO-K Tray, L-4-S	2 11/16	68 mm	±3/16	±5 mm	2 8/16 to 2 14/16	63 - 73 mm	2 7/8	73 mm	±3/16	±5 mm	2 11/16 to 3 1/16	68-78 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE L5S	263485	Perf. STO-K Tray, L-5-S	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 3/4	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65-75 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE U1T	263571	Perf. STO-K Tray, U-1-T	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2 1/8	54 mm	±3/16	±5 mm	1 15/16 to 2 5/16	49-59 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE U2T	263572	Perf. STO-K Tray, U-2-T	2 5/8	67 mm	±3/16	±5 mm	2 7/16 to 2 13/16	62 - 72 mm	2 1/4	57 mm	±3/16	±5 mm	2 1/16 to 2 7/16	52-62 mm	7/16	11 mm	±3/16	±5 mm	4/16 to 10/16	6-16 mm
COE U3T	263573	Perf. STO-K Tray, U-3-T	3 1/8	79 mm	±3/16	±5 mm	2 15/16 to 3 5/16	74 - 84 mm	2 11/16	68 mm	±3/16	±5 mm	2 8/16 to 2 14/16	63-73 mm	1/2	13 mm	±3/16	±5 mm	5/16 to 11/16	8-18 mm
COE L1T	263581	Perf. STO-K Tray, L-1-T	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	56 - 66 mm	2 9/16	65 mm	±3/16	±5 mm	2 6/16 to 2 12/16	60-70 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE L2T	263582	Perf. STO-K Tray, L-2-T	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2 5/16	59 mm	±3/16	±5 mm	2 2/16 to 2 8/16	54-64 mm	1/2	13 mm	±3/16	±5 mm	5/16 to 11/16	8-18 mm
COE L3T	263583	Perf. STO-K Tray, L-3-T	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57 - 67 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	1/2	13 mm	±3/16	±5 mm	5/16 to 11/16	8-18 mm
COE L4T	263584	Perf. STO-K Tray, L-4-T	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	3/8	10 mm	±3/16	±5 mm	3/16 to 9/16	5-15 mm
COE U10	263671	Perf. STO-K Tray, U-1-0	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59 - 69 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	9/16	14 mm	±3/16	±5 mm	6/16 to 12/16	9-19 mm
COE U20	263672	Perf. STO-K Tray, U-2-0	2 5/8	67 mm	±3/16	±5 mm	2 7/16 to 2 13/16	62 - 72 mm	2 1/4	57 mm	±3/16	±5 mm	2 1/16 to 2 7/16	52-62 mm	5/8	16 mm	±3/16	±5 mm	7/16 to 13/16	11-21 mm
COE U30	263673	Perf. STO-K Tray, U-3-0	2 12/16	70 mm	±3/16	±5 mm	2 9/16 to 2 15/16	65 - 75 mm	2 3/8	60 mm	±3/16	±5 mm	2 3/16 to 2 9/16	55-65 mm	1/2	13 mm	±3/16	±5 mm	5/16 to 11/16	8-18 mm
COE U40	263674	Perf. STO-K Tray, U-4-0	3	76 mm	±3/16	±5 mm	2 13/16 to 3 3/16	71 - 81 mm	2 1/2	64 mm	±3/16	±5 mm	2 5/16 to 2 11/16	59-69 mm	11/16	17 mm	±3/16	±5 mm	8/16 to 14/16	12-22 mm
COE L10	263681	Perf. STO-K Tray, L-1-0	2 7/16	62 mm	±3/16	±5 mm	2 4/16 to 2 10/16	57 - 67 mm	2	51 mm	±3/16	±5 mm	1 13/16 to 2 3/16	46-56 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE L20	263682	Perf. STO-K Tray, L-2-0	2 11/16	68 mm	±3/16	±5 mm	2 8/16 to 2 14/16	63 - 73 mm	2 1/8	54 mm	±3/16	±5 mm	1 15/16 to 2 5/16	49-59 mm	5/16	8 mm	±3/16	±5 mm	2/16 to 8/16	3-13 mm
COE L30	263683	Perf. STO-K Tray, L-3-0	2 15/16	75 mm	±3/16	±5 mm	2 12/16 to 3 2/16	7												