



GC AMERICA INC.  
 3737 W 127<sup>th</sup> STREET  
 ALSIP, ILLINOIS 60803  
 TEL (800) 323-3386  
 FAX (708) 897-4062  
 www.gcamerica.com

**Contact:**

Toni M. Cusimano  
 Assoc. Laboratory Product Mgr. & CE Mgr.  
 GC America Inc.  
 800-323-7063 ext. 3006

**FOR IMMEDIATE RELEASE:**

Alsip, IL – October 2016 - GC America Inc. introduces GC LiSi™ PressVest, a phosphate bonded investment for press ceramic techniques offering you an easier, faster and safer process. Hydrofluoric acid treatment and alumina blasting are no longer required for the removal of the reaction layer. GC LiSi PressVest provides a stable handling under various temperature conditions and is so robust it provides a strong investment ring so you can effectively use wider volume, but also easy to divest. GC LiSi PressVest provides a comfortable and secure investing with a high flow and a longer working time.

**Product will be available to ship on November 14, 2016**  
**Canadian launch to follow**

- Easy removable of the reaction layer
  - Reduction of operating time
  - Reproducible fit due to no use of hydrofluoric acid
- Easy divesting
  - Reduction of operating time
- Super high fluidity and long working time
  - Accurate and detailed reproduction
- Precise expansion control
  - Optimal fit
- Optimized for quick heating technique
- Adapted to multiple press ceramic technique

SKU#	Product Description
901424	GC LiSi PressVest Powder 100 g x 60
901425	GC LiSi PressVest Liquid 900mL
901426	GC LiSi PressVest SR (Surface Refining) Liquid 100mL
901427	GC LiSi PressVest Intro Kit (Kit contents: (6) LiSi PressVest Powder 100g, (1) LiSi PressVest Liquid, 135mL, (1) LiSi PressVest SR (Surface Refining) Liquid 30mL, (1) Spray Nozzle) <b>PLEASE NOTE: KIT CONTENTS ARE NOT AVAILABLE FOR INDIVIDUAL PURCHASES</b>

###

