



# MultiPressVest from GC.



## Carbon-Free Phosphate Bonded Investment For Pressable Ceramics

GC America has specifically developed a new investment material adapted to multiple press ceramic techniques offered within the GC Initial PC pressable system and the GC Initial IQ – One Body, Press-over-Metal/Zircon system.

### Characteristics & benefits

- Fine-grained phosphate bonded investment produces smooth and detailed pressings
- Creamy consistency ensures easy and accurate investing
- High thermal stability makes it suitable for the quick heating as well as for the conventional heat-up technique

### Specially developed for multiple press ceramic techniques like

- All ceramic press techniques with Initial PC: crowns, inlays and facings
- Press ceramic over metal structures with GC Initial IQ – One Body, Press-over-Metal.
- Press ceramic over zirconium structures with GC Initial IQ – One Body, Press-over-Zircon.

Art code	Description
800242	GC MultiPressVest, Intro pack (powder 6 x 100 g bags, 135 ml liquid)
800243	GC MultiPressVest, 6 kg powder (60 x 100 g bags)
800244	GC MultiPressVest, 900 ml liquid



GC AMERICA INC.  
 3737 WEST 127TH STREET  
 ALSIP, ILLINOIS 60803  
 TEL (708) 597-0900  
 FAX (708) 371-5103