

GC AMERICA'S CANADIAN JANUARY-MARCH 2018 PROMOTIONS

RESTORATIVES

Glass Hybrid

- Buy 4 EQUIA® Forte Products, **Get 1 FREE!***

Glass Ionomer

- Buy 4 Glass Ionomer Restorative Products, **Get 1 FREE!***



Composite/Bonding

- Buy 4 G-aenial Sculpt® Refills or G-Premio BOND™ Products, **Get 1 FREE!***
- Buy 4 G-aenial™ Flo or Universal Flo Refills, **Get 1 FREE!***



In-Office Milling**

- Buy 3 CERASMART™ Packages, **Get 1 G-CEM LinkForce™ Starter Kit FREE!***
- Buy 2 CERASMART™ Packages, **Get 1 G-CEM LinkForce™ Refill FREE!***
- Buy 5 CERASMART™ Packages and 1 Complete OPTIGLAZE™ color Kit, **Get 1 The Light 405 FREE!***
- Buy 5 CERASMART™ Packages, **Get 1 FREE!***
- Buy 6 CERASMART™ Packages, **Get 1 CERAMIC PRIMER II & 1 OPTIGLAZE™ FREE!***
- Buy 2 CERASMART™ Packages, **Get 1 GRADIA® DiaPolisher FREE!***

EQUIPMENT



- Buy 6 Capsule Intro Kits or Refills, **Get 1 CAPSULE MIXER CM-II™ FREE!***
- Buy \$1,200 of Restorative and Cement Products, **Get 1 The Light or The Light 405 FREE!***
Excludes GC Fuji® powder/liquid, GC Fuji ORTHO™ products, G-COAT PLUS™, GC ORTHO GEL CONDITIONER and UniFil® Bond products.
- Buy \$2,250 of GC Products, **Get 1 Renfert® SYMPRESS Mixing Machine FREE!***
Contact GC America Customer Service for a list of qualifying products.

CEMENTS

Resin

- Buy 3 G-CEM LinkForce™ Refills, **Get 1 FREE!***
- Buy 2 G-CEM LinkForce™ Refills, **Get 1 CERASMART™ Package FREE!***
- Buy 3 G-CEM LinkAce™ Refills, **Get 1 FREE!***

Resin Modified Glass Ionomer

- Buy 4 GC Fuji CEM™ 2 or GC Fuji PLUS® Products, **Get 1 FREE!***



PREVENTIVE



- Buy 3 MI Paste™ or MI Paste Plus™ 10-Tube Packages, **Get 1 FREE!***
- Buy 3 MI Varnish™ Packages, **Get 1 FREE!***

IMPRESSION MATERIALS

- Buy 4 EXAFAST™, EXAMIX™ or EXAJET™ Products, **Get 1 FREE!***
Excludes 32-packs
- Buy 4 EXA'lence™ 370 Products, **Get 1 FREE!***
- Buy 4 Tray Products, **Get 1 FREE!***



TEMPORIZATION



- Buy 4 TEMPSMART® Products, **Get 1 FREE!***

GC America Special Offers Redemption Information:

You can redeem your free goods by phone or fax. An authorized dealer invoice can be sent to GC America Inc. Attn: Customer Care Group, 3737 West 127th Street, Alsip, IL 60803 or faxed to: 800.423.2963.

PHONE _____

DOCTOR NAME _____

ADDRESS _____


CITY _____ STATE _____ ZIP _____



STRENGTH & ESTHETICS IN ONE SYSTEM

G-CEM LinkForce™ Dual-Cure Adhesive Resin Cement

- G-CEM LinkForce is a strong, esthetic, universal solution for all adhesive cementations
- G-CEM LinkForce is the ideal cement to use with our **NEW GC Initial™ LiSi Press Lithium Disilicate Ingots**
- Four shades available in both cement and try-in paste
- Excellent color stability and easy clean-up



Source: GCC R&D
†Not a registered trademark of GC America Inc.



Buy 3 G-CEM LinkForce Refills,
Get 1 FREE!*



Use G-CEM LinkForce to cement your GC Initial LiSi Press Lithium Disilicate crowns!

GC FujiCEM™ 2 Resin Modified Glass Ionomer Cement

Buy 4,
Get 1 FREE!*



- Strong Retention: Flex Fuse Technology
- Assured Sealing Ability: Hyper-mineralized zone at tooth interface
- Stress-free Seating: Moisture tolerant, unmatched film thickness with a creamy consistency
- No Post-operative Sensitivity: Biocompatible to tooth structure and surrounding soft tissue
- Minimized Risk of Secondary Caries: Sustained fluoride release

G-CEM LinkAce™ Self-Adhesive Resin Cement

Buy 3 Refills,
Get 1 FREE!*



- Optimal Self-cure Mode: Innovative initiator system (chemical)
- High Bond Durability to Zirconia in ONE Step: Proprietary phosphate monomers
- Exceptional Color Stability: Very low water sorption and HEMA-free
- Unsurpassed Wear Resistance: Small homogeneously distributed particles
- Easy Excess Cement Removal: Only 1-2 second tack cure