

CERAMICS



- Buy 3 GC Initial® LiSi Press Ingot Sleeves, **Get 3 GC Initial® LiSi 20g Refills FREE!***
- Buy 3 GC Initial® LiSi Press Ingot Sleeves, **Get 3 GC Initial® IQ Lustre Pastes Refills OR GC Initial Spectrum Stains FREE!***
- Buy 5 GC Initial® LiSi Press Ingot Sleeves, **Get 1 GC LiSi PressVest Powder & Liquid FREE!***
- Buy 2 GC Initial® LiSi Press Ingot Sleeves, **Get 1 GC LiSi PressVest Intro Kit FREE!***
- Buy 2 GC Initial® Sets, **Get 1 FREE!***
- Buy 2 GC Initial® Zirconia Disks, **Get 1 FREE!***
- Buy 1 GC Initial® Set, **Get 1 GC Initial® Brush Set FREE!***
- Buy 6 GC Initial® Refills, **Get 1 GC Initial® Brush FREE!***
- Buy 3 GC Initial® Refills, **Get 1 of any GC Initial® Brush Tips FREE!***

DENTURE BASE RESINS

- Buy 1 NATURE-CRYL® 25lb Powder, **Get 4 NATURE-CRYL® 32oz Liquids FREE!***
- Buy 1 NATURE-CRYL® 5lb Powder, **Get 1 NATURE-CRYL® 32oz Liquid FREE!***



CAD/CAM



- Buy 2 GC Initial® LRF Packages, **Get 1 GC Initial® LRF Glaze Paste & Liquid FREE!***
- Buy 3 GC Initial® LRF Packages, **Get 1 FREE!***
- Buy 3 CERASMART® Packages, **Get 1 FREE!***
- Buy 6 CERASMART® Packages, **Get 1 G-CEM LinkForce® System Kit FREE!***
- Buy 6 GC Initial® LRF Packages, **Get 1 G-CEM LinkForce® System Kit FREE!***
- Buy 5 CERASMART® Packages & 1 OPTIGLAZE™ color Complete Set, **Get 1 The Light 405 FREE!***

INVESTMENTS

- Buy any GC Investment Powder Product (GC Fujivest® II, GC Fujivest® Premium, GC Fujivest® Platinum II, GC MultiPressVest®, GC LiSi PressVest), **Get 1 Investment 900mL Liquid FREE!***



Weather temperature must remain over 32°F in Chicago land area

COMPOSITES



- Buy 1 GRADIA® gum shades Starter Kit, **Get 1 CERAMIC PRIMER II, METALPRIMER II or COMPOSITE PRIMER FREE!***
- Buy 3 OPTIGLAZE™ color Refills, **Get 1 FREE!***

GYPSUMS & MODELING

- Buy 4 GC FUJIROCK® Packages, **Get 1 FREE!***

**GC America Special Offers Redemption Information:**

You can redeem your free goods by phone or fax. An authorized dealer invoice can be sent to GC America Inc. Attn: Customer Care Group, 3737 West 127th Street, Alsip, IL 60803 or faxed to: 800.423.2963.

PHONE _____

DOCTOR/OFFICE NAME _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____



Natural beauty restored.



Courtesy of Al Hodges, CDT

Imagine a pressable ceramic that is exceptionally strong, extremely durable, with amazing aesthetics and saves your lab significant time!

- Excellent flexural strength (>500 MPa*) with HDM Technology
- Save up to 20 minutes per restoration
- Amazing aesthetics
- Rich, warm, and bright colors with excellent fluorescence
- Virtually no reaction layer when divested

*Data on file.



Lithium Disilicate Redefined

