

CERAMICS



- Buy 4 of ANY GC Initial® Refills, **Get 1 FREE!*** same type
- Buy 5 GC Initial® LiSi Press Ingot Sleeves, **Get 1 GC LiSi Press Vest Powder & Liquid FREE!***
- Buy 2 GC Initial® Sets, **Get 1 FREE!***
- Buy 2 GC Initial® Zirconia Disks, **Get 1 FREE!***
- Buy 1 GC Initial® Set, **Get 1 GC Initial® Brush Set FREE!***
- Buy 6 GC Initial® Refills, **Get 1 GC Initial® Brush FREE!***
- Buy 3 GC Initial® Refills, **Get 1 of any GC Initial® Brush Tips FREE!***

COMPOSITES



- Buy 2 GRADIA® PLUS Kits, **Get 1 FREE!***
- Buy 3 GRADIA® PLUS Refills, **Get 1 FREE!*** includes primer
- Buy 1 Labolight DUO, **Get 1 GRADIA® PLUS Kit FREE!***
- Buy 1 GRADIA® PLUS Kit, **Get 1 GRADIA® PLUS Instrument Kit FREE!***
- Buy 3 OPTIGLAZE™ color Refills, **Get 1 FREE!***

DENTURE BASE RESINS

- Buy 1 NATURE-CRYL® 25lb Powder, **Get 4 NATURE-CRYL® 32oz Liquids FREE!*** As weather permits
- Buy 1 NATURE-CRYL® 5lb Powder, **Get 1 NATURE-CRYL® 32oz Liquid FREE!*** As weather permits



CAD/CAM



- Buy 2 GC Initial® LRF Packages, **Get 1 GC Initial® LRF Glaze Paste & Liquid FREE!***
- Buy 3 GC Initial® LRF Packages, **Get 1 FREE!***
- Buy 3 CERASMART® 270 Packages, **Get 1 FREE!***

3D PRINTING

- Buy 4 GC Temp PRINT™ 500g Bottles, **Get 1 OPTIGLAZE™ color Kit FREE!***



GYPSUMS & MODELING

- Buy 3 GC FUJIROCK® Packages, **Get 1 FREE!***



GC America Special Offers Redemption Information:

You can redeem your free goods by e-mail, mail, phone or fax. An authorized dealer invoice can be sent to:

- promotions.gcamerica@gc.dental
- GC America Inc. Attn: Customer Care Group, 3737 West 127th Street, Alsip, IL 60803
- Call 800.323.7063
- Faxed to 800.423.2963

PHONE _____

DOCTOR/OFFICE NAME _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____



Since 1921
100 years of Quality in Dental



BUILDING BRIDGES FOR THE FUTURE

GC Temp PRINT™

3D Printable Light-Cured Composite for Temporary Crown and Bridge

Unique Filler Technology

- Homogeneous dispersion of surface coated silica fillers lend to exceptional wear resistance
- Temporary restorations with high resistance to abrasion in the oral environment
- Ideal for long-term provisional crown and bridges

Dynamic Controlled Rheology (DCR)

- Anti-sedimentation additives ensure the highest precision and material stability
- Anti-sedimentation additives form a shell around the pigments and silica-fillers, preventing precipitation
- Anti-sedimentation additives allow for simple handshaking to sufficiently obtain a homogeneous dispersion, even after long-term storage



Buy 4 GC Temp PRINT™ 500g Bottles, **Get 1 OPTIGLAZE™ color Kit FREE!***

800.323.7063 • www.gcamerica.com • www.gcamerica.com/training • customerservice.gcamerica@gc.dental • © 2021 GC America Inc.

*Promo code 2519. Terms and Conditions of ALL Promotions: Limit three redemptions per customer without prior written authorization. All offers expire 12/31/2021. Effective Dates: 7/1/2021 - 12/31/2021. Promotions cannot be combined with any other promotional offers. No substitutions allowed. Offers are valid in the US and Canada only. Free goods must be redeemed within six weeks of dealer invoice date. Free goods will be of equal or lesser value to the lowest selling unit price. Please allow 4-8 weeks for free goods delivery. All free goods will be fulfilled by GC America Inc. except where noted. Please indicate requested free goods on dealer invoice. Only free good requests that come directly from the Dental Office/Dentist and/or Staff will be honored. Qualifying purchases must appear on a single dealer invoice.